


LePuffer™

Single-use Balloon Dilatation Catheter
Progressive Dilation Balloons




Leo+Med
乐奥医疗 为生命健康服务

GI Balloon Dilatation Catheter Series 

LePuffer™ Single-use Balloon Dilatation Catheter Progressive Dilation Balloons

Used to dilate strictures of the biliary duct and the alimentary tract.



• Excellent Deliverability

- Soft and tapered tip allows less invasive advancement across the stricture.
- Low-profile design reduces resistance and increases crossability and trackability.



• Operational Efficiency

- Facilitate placement within tight strictures and tortuous anatomy.
- Optimized balloon deflation time.



• Optimized Balloon Design

- The balloon component is designed to offer stability during stricture dilation.
- Semi-compliant material engineered to adapt to biliary duct.

Specifications

	REF	Balloon			Catherer			Minimum Working Channel (mm)
		Length (mm)	Diam. (mm)	Nominal Pressure (atm)	Length (cm)	Diam. (mm)	Structure	
Pyloric and Biliary	KZ0630E0	30	6-7-8	3-4.5-6	180	2	OTW	2.8
	KZ0830E0	30	8-9-10	3-4.5-6	180	2	OTW	2.8
	KZ1030E0	30	10-11-12	3-4.5-6	180	2	OTW	2.8
	KZ1230E0	30	12-13.5-15	3-4.5-6	180	2	OTW	2.8
	KZ1530E0	30	15-16.5-18	3-4.5-6	180	2	OTW	2.8
Esophageal, Achalasia, Pyloric and Biliary	KZ0655G0	55	6-7-8	3-4.5-6	180	2.5	OTW	2.8
	KZ0855G0	55	8-9-10	3-4.5-6	180	2.5	OTW	2.8
	KZ1055G0	55	10-11-12	3-4.5-6	180	2.5	OTW	2.8
	KZ1255G0	55	12-13.5-15	3-4.5-6	180	2.5	OTW	2.8
	KZ1555G0	55	15-16.5-18	3-4.5-6	180	2.5	OTW	2.8
	KZ1855G0	55	18-19-20	3-4.5-5	180	2.5	OTW	2.8
	KZ2055G0	55	20-21-22	3-4.5-6	180	2.5	OTW	2.8
	KZ1280F0	80	12-13.5-15	3-4.5-6	180	2.5	OTW	2.8
	KZ1580F0	80	15-16.5-18	3-4.5-6	180	2.5	OTW	2.8
	KZ1880F0	80	18-19-20	3-4.5-6	180	2.5	OTW	2.8
Colonic	KZ0655H0	55	6-7-8	3-4.5-6	240	2.5	OTW	2.8
	KZ0855H0	55	8-9-10	3-4.5-6	240	2.5	OTW	2.8
	KZ1055H0	55	10-11-12	3-4.5-6	240	2.5	OTW	2.8
	KZ1255H0	55	12-13.5-15	3-4.5-6	240	2.5	OTW	2.8
	KZ1555H0	55	15-16.5-18	3-4.5-6	240	2.5	OTW	2.8
	KZ1855H0	55	18-19-20	3-4.5-6	240	2.5	OTW	2.8
	KZ2055H0	55	20-21-22	3-4.5-5	240	2.5	OTW	2.8

▲ RX Structure Catheter is also available. Please consult with LEOMED representative for details.